

PDQ

Interference
Free

A_{1c}

in 2
minutes

the
definitive
technology
for
GHb/A_{1c}
testing



- ◆ Compact HPLC system
- ◆ Proven accurate boronate affinity technology
- ◆ Quick: answers in 2 minutes
- ◆ Interference-free results
- ◆ Flexible for low or moderate volume testing
- ◆ Simple, menu driven software
- ◆ Lightweight, minimum footprint
- ◆ Low initial capital investment
- ◆ Minimal maintenance
- ◆ Ideal for central or satellite laboratories and clinics



- ◆ The **PDQ PLUS** AutoSampler (shown at right) provides primary tube sampling with on-line barcode reading.
- ◆ Primus Reagents have 2-year room temperature stability.
- ◆ LIS Interfaceable.
- ◆ Disk storage of test results.
- ◆ Built-in diagnostic monitoring system.
- ◆ Complete, individual test reports and test summary from each batch.



- ◆ The Primus **PDQ** is also available as a stand-alone unit that allows your laboratory to process individual samples.

<p>SOFTWARE PACKAGE: Patient Identification – 9 digit keypad entry or barcode Patient data storage – 200 results User selectable report format Automatic or manual calibration Diagnostic alerts to assist operator Levey-Jennings plot of QC data</p>	<p>AMBIENT LIMITS: Temperature – 15 - 32°C (60 to 90°F) Relative Humidity – 10 - 90% non-condensing</p> <p>DIMENSIONS: Height – 34 cm (14 in) Width – 45 cm (18 in) Depth – 42 cm (17 in) Weight – 19 Kg (42 lb)</p>
---	--



7235 Central • Kansas City, MO 64114 • www.primusdiagnostics.com

816-523-7491 • 800-377-4752 • 816-361-1974 FAX